

## Anti-Psmd11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00231

Item Number: CM0015911



### Introduction

Rabbit polyclonal antibody targeting PSMD11, a 26S proteasome subunit crucial for assembly. Suitable for WB, IF/ICC, and ELISA applications in human, mouse, and rat samples.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | PSMD11 Rabbit pAb                                                                                                    |
| Remarks/Alias    | S9; Rpn6; p44.5; PSMD11                                                                                              |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 5717                                                                                                                 |
| Gene ID          | 5717                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 253-422 of human PSMD11 (NP_002806.2). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, IF/ICC, ELISA                                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | O00231                                                                                                               |
| Protein Weight   | 47kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |