


Anti-Psmd13 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Unm6

Item Number: CM0030052



Introduction

High-quality anti-PSMD13 rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat PSMD13, a 26S proteasome regulatory subunit crucial for protein degradation. Ideal for proteasome and ubiquitin-proteasome pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PSMD13 Rabbit pAb
Remarks/Alias	S11; Rpn9; p40.5; HSPC027; PSMD13
Species	Human
Gene ID (Human)	5719
Gene ID	5719
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human PSMD13 (NP_002808.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UNM6
Protein Weight	43kDa
Shipping	Ice bag