

Anti-Psmd14 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O00487

Item Number: CM0014801



Introduction

Rabbit polyclonal antibody against human PSMD14, a 26S proteasome subunit. Suitable for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Critical for ubiquitin-dependent protein degradation, DNA repair, and autophagy.

[Learn More](#)

Parameter	Value
Product Name	PSMD14 Rabbit pAb
Remarks/Alias	PAD1; POH1; RPN11; PSMD14
Species	Human
Gene ID (human)	10213
Gene ID	10213
Immunogen	Recombinant fusion protein containing amino acids 160-300 of human PSMD14 (NP_005796.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00487
Protein Weight	35kDa
Shipping	Ice bag