

Anti-Psmd2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13200

Item Number: CM0019119



Introduction

High-quality rabbit polyclonal antibody targeting human PSMD2 (26S proteasome regulatory subunit S2, p97), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for studying ubiquitin-mediated proteolysis, cell cycle, and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	PSMD2 Rabbit pAb
Remarks/Alias	S2, P97, RPN1, TRAP2, PSMD2
Species	Human
GeneID (Human)	5708
GeneID	5708
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PSMD2 (NP_002799.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13200
Protein Weight	100kDa
Shipping	Ice bag