


Anti-Psmd3 Rabbit Polyclonal Antibody For Wb, Elisa - O43242

Item Number: CM0028681



Introduction

Rabbit polyclonal antibody against human PSMD3, a 26S proteasome regulatory subunit involved in ubiquitinated protein degradation. Validated for WB and ELISA, cross-reacts with mouse and rat, suitable for proteasome and protein homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	PSMD3 Rabbit pAb
Remarks/Alias	S3, P58, RPN3, TSTA2, PSMD3
Species	Human
GeneID (Human)	5709
GeneID	5709
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43242
Protein Weight	61kDa
Shipping	Ice bag