


Anti-Psmd4 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P55036

Item Number: CM0010555



Introduction

High-specificity rabbit monoclonal antibody against PSMD4 (26S proteasome subunit S5A/Rpn10), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Ideal for ubiquitin-proteasome pathway research and proteasome dysfunction studies.

[Learn More](#)

Parameter	Value
Product Name	PSMD4 Rabbit mAb
Remarks/Alias	AF, ASF, S5A, AF-1, MCB1, Rpn10, pUB-R5, PSMD4
Species	Human
GeneID (Human)	5710
GeneID	5710
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human PSMD4 (P55036).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55036
Protein Weight	41kDa
Shipping	Ice bag