


Anti-Psmd4 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P55036

Item Number: CM0015845



Introduction

PSMD4 Rabbit pAb is a polyclonal antibody from rabbit, validated for WB, IF/ICC, IP, and ELISA. Reacts with human and rat PSMD4, a 26S proteasome ubiquitin receptor subunit, supporting ubiquitin-proteasome research.

[Learn More](#)

Parameter	Value
Product Name	PSMD4 Rabbit pAb
Remarks/Alias	AF; ASF; S5A; AF-1; MCB1; Rpn10; pUB-R5; PSMD4
Species	Human
GeneID (Human)	5710
GeneID	5710
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	P55036
Protein Weight	41kDa
Shipping	Ice bag