


Anti-Psmd6 Polyclonal Antibody For Wb, Elisa - Q15008

Item Number: CM0022208



Introduction

Rabbit polyclonal antibody against human PSMD6 (26S proteasome regulatory subunit S10/Rpn7). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for proteasome and protein homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	PSMD6 Rabbit pAb
Remarks/Alias	S10; Rpn7; p42A; p44S10; SGA-113M; PSMD6
Species	Human
Gene ID (Human)	9861
Gene ID	9861
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-307 of human PSMD6 (NP_055629.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15008
Protein Weight	46kDa
Shipping	Ice bag