

Anti-Psmd7 Rabbit Polyclonal Antibody - P51665

Item Number: CM0026777



Introduction

Rabbit polyclonal antibody against human PSMD7 (26S proteasome regulatory subunit) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse. Suitable for proteasome and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	PSMD7 Rabbit pAb
Remarks/Alias	P40; S12; Rpn8; MOV34; PSMD7
Species	Human
GeneID (Human)	5713
GeneID	5713
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-324 of human PSMD7 (NP_002802.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51665
Protein Weight	37kDa
Shipping	Ice bag