

Anti-Psmd9 Polyclonal Antibody For Wb, Elisa - 000233

Item Number: CM0030029

Introduction

Recombinant-based rabbit polyclonal antibody targeting PSMD9 (p27/Rpn4). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Investigate proteasome assembly and chaperone function.

[Learn More](#)

Parameter	Value
Product Name	PSMD9 Rabbit pAb
Alias	p27, Rpn4, PSMD9
Species (Target)	Human
Gene ID	5715
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-223 of human PSMD9 (NP_002804.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00233
Protein Weight	25kDa
Shipping	Ice bag