


Anti-Psme1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q06323

Item Number: CM0025302



Introduction

Rabbit polyclonal antibody against PSME1 (PA28alpha), a proteasome activator involved in antigen processing. Validated for WB, IF/ICC, ELISA in human and mouse samples.

Immunogen: recombinant fusion protein of human PSME1 (aa 1-249).

[Learn More](#)

Parameter	Value
Product Name	PSME1 Rabbit pAb
Remarks/Alias	PA28A; IFI5111; REGalpha; PA28alpha; HEL-S-129m; PSME1
Species	Human
GeneID (Human)	5720
GeneID	5720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-249 of human PSME1 (NP_006254.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06323
Protein Weight	29kDa
Shipping	Ice bag