


## Anti-Psme1/Pa28A Rabbit Monoclonal Antibody For Wb, Elisa - Q06323

Item Number: CM0010822



### Introduction

Rabbit monoclonal antibody targeting human PSME1 (PA28 $\alpha$ ). Suitable for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying immunoproteasome function and antigen processing.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | PSME1 Rabbit mAb                                                                                                 |
| Remarks/Alias    | PA28A; IFI5111; REGalpha; PA28alpha; HEL-S-129m; PSME1                                                           |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 5720                                                                                                             |
| GeneID           | 5720                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human PSME1 (NP_006254.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q06323                                                                                                           |
| Protein Weight   | 29kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |