


Anti-Psme2 Polyclonal Antibody (Ko Validated) For Wb And Elisa - Q9UI46

Item Number: CM0023186



Introduction

Rabbit polyclonal antibody against human PSME2 (PA28 β), validated for Western blot and ELISA. Reacts with human and mouse; targets the immunoproteasome activator involved in antigen processing.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PSME2 Rabbit pAb
Remarks/Alias	PA28B; REGbeta; PA28beta; E2
Species	Human
Gene ID (Human)	5721
Gene ID	5721
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human PSME2 (NP_002809.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UL46
Protein Weight	27kDa
Shipping	Ice bag