


Anti-Psme2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9UL46

Item Number: CM0024839



Introduction

KO validated rabbit polyclonal antibody against human PSME2 (PA28 β) for western blotting and ELISA applications; cross-reacts with mouse. Suitable for immunoproteasome and antigen processing research.

[Learn More](#)

Product Name	[KO Validated] PSME2 Rabbit pAb
Remarks/Alias	PA28B; REGbeta; PA28beta; E2
Species	Human
GeneID (human)	5721
GeneID	5721
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human PSME2 (NP_002809.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UL46
Protein Weight	27kDa
Shipping	Ice bag