

Anti-Psme3 Polyclonal Antibody For Wb, Ihc-P - P61289

Item Number: CM0014182



Introduction

High-quality rabbit polyclonal antibody against human PSME3 (PA28gamma, REG-gamma), validated for WB, IHC-P, IF/ICC, ELISA. Detects endogenous PSME3 in human, mouse, rat. Ideal for proteasome and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	PSME3 Rabbit pAb
Remarks/Alias	Ki; PA28G; HEL-S-283; PA28gamma; REG-GAMMA; PA28-gamma; PSME3
Species	Human
Gene ID (Human)	10197
Gene ID	10197
Immunogen	Synthetic peptide corresponding to amino acids 200-254 of human PSME3 (NP_005780.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61289
Protein Weight	30kDa
Shipping	Ice bag