

Anti-Psme3 Rabbit Polyclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Elisa - P61289

Item Number: CM0018913



Introduction

This [KO Validated] PSME3 rabbit polyclonal antibody detects human, mouse, and rat PSME3, validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for proteasome and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PSME3 Rabbit pAb
Remarks/Alias	Ki; PA28G; HEL-S-283; PA28gamma; REG-GAMMA; PA28-gamma; E3
Species	Human
GeneID (Human)	10197
GeneID	10197
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-254 of human PSME3 (NP_005780.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61289
Protein Weight	30kDa
Shipping	Ice bag