


Anti-Psme4 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14997

Item Number: CM0016748



Introduction

Rabbit polyclonal antibody against PSME4 (PA200), validated in WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Ideal for studying proteasome-mediated histone degradation, spermatogenesis, and DNA damage response.

[Learn More](#)

Parameter	Specification
Product Name	PSME4 Rabbit pAb
Remarks/Alias	PA200; PSME4
Species	Human
Gene ID (Human)	23198
Gene ID	23198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1634-1843 of human PSME4 (NP_055429.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14997
Protein Weight	211kDa
Shipping	Ice bag