


Anti-Psmf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q92530

Item Number: CM0026747



Introduction

PSMF1 Rabbit Polyclonal Antibody designed for Western blot (WB), IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat PSMF1. Proteasome inhibitor PI31 subunit, key in proteasome function control.

[Learn More](#)

Parameter	Value
Product Name	PSMF1 Rabbit pAb
Remarks/Alias	PI31; PSMF1
Species	Human
GeneID (Human)	9491
GeneID	9491
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-271 of human PSMF1 (NP_006805.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92530
Protein Weight	30kDa
Shipping	Ice bag