


Anti-Psmg1 Polyclonal Antibody For Wb, Ihc-P, Elisa - 095456

Item Number: CM0017487



Introduction

Rabbit polyclonal antibody against human Proteasome assembly chaperone 1 (PSMG1); validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for proteasome assembly and Down syndrome research.

[Learn More](#)

Parameter	Value
Product Name	PSMG1 Rabbit pAb
Remarks/Alias	PAC1, DSCR2, PAC-1, C21LRP, LRPC21, PSMG1
Species	Human
GeneID	8624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-288 of human PSMG1 (NP_003711.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	095456
Protein Weight	33kDa
Shipping	Ice bag