


Anti-Pspc1 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q8Wxf1

Item Number: CM0029949



Introduction

Rabbit polyclonal antibody targeting human PSPC1, a paraspeckle component and RNA-binding protein. Validated for WB, IHC-P, IP, and ELISA. Cross-reacts with mouse and rat. Ideal for paraspeckle and transcription research.

[Learn More](#)

Parameter	Value
Product Name	PSPC1 Rabbit pAb
Remarks/Alias	PSP1; PSPC1
Species	Human
Gene ID (Human)	55269
Gene ID	55269
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human PSPC1 (NP_001035879.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXF1
Protein Weight	59kDa
Shipping	Ice bag