

Anti-Psph Rabbit Polyclonal Antibody For Wb And Elisa - P78330

Item Number: CM0030463



Introduction

Rabbit polyclonal antibody against human phosphoserine phosphatase (PSPH) validated for Western blot and ELISA, reactive with human, mouse, rat. For serine biosynthesis research. P78330.

[Learn More](#)

Parameter	Value
Product Name	PSPH Rabbit pAb
Remarks/Alias	PSP; PSPHD; PSPH
Species	Human
Gene ID (Human)	5723
Gene ID	5723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human PSPH (NP_004568.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78330
Protein Weight	25kDa
Shipping	Ice bag