

Anti-Psrc1 Polyclonal Antibody For Wb, Elisa - Q6Pgn9

Item Number: CM0021108

Introduction

High-quality rabbit polyclonal antibody targeting human PSRC1 (DDA3), validated for WB and ELISA. Cross-reacts with mouse. Ideal for mitotic spindle and chromosome segregation research.

[Learn More](#)

Parameter	Value
Product Name	PSRC1 Rabbit pAb
Remarks/Alias	DDA3; FP3214; PSRC1
Species	Human
Gene ID (Human)	84722
Gene ID	84722
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human PSRC1 (NP_001005290.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6PGN9
Protein Weight	39kDa
Shipping	Ice bag