


Anti-Pstpip1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O43586

Item Number: CM0026494



Introduction

High-quality rabbit polyclonal antibody against human PSTPIP1, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying actin cytoskeleton regulation, T-cell activation, and inflammasome biology.

[Learn More](#)

Parameter	Value
Product Name	PSTPIP1 Rabbit pAb
Remarks/Alias	H-PIP; PAPAS; CD2BP1; PSTPIP; CD2BP1L; CD2BP1S; PSTPIP1
Species	Human
GeneID (Human)	9051
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43586
Protein Weight	48kDa
Shipping	Ice bag