

Anti-Ptbp1 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P26599

Item Number: CM0025402



Introduction

PTBP1 Rabbit pAb (P26599) detects human, mouse, rat PTBP1 by WB, IF/ICC, IHC-P, ELISA. This polyclonal antibody targets the polypyrimidine tract-binding protein 1, a key regulator of alternative splicing, essential for RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	PTBP1 Rabbit pAb
Remarks/Alias	PTB; PTB2; PTB3; PTB4; pPTB; HNRPI; PTB-1; PTB-T; HNRNPI; HNRNP-I; PTBP1
Species	Human
Gene ID (Human)	5725
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human PTBP1 (NP_114368.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26599
Protein Weight	60 kDa
Shipping	Ice bag