


Anti-Ptbp2 Rabbit Monoclonal Antibody For Wb, Ip, If And Ihc - Q9Uka9

Item Number: CM0005638



Introduction

Rabbit monoclonal antibody against human PTBP2, validated for Western blot, immunoprecipitation, immunofluorescence, and immunohistochemistry; reacts with human, mouse, rat PTBP2. Ideal for splicing and translation research.

[Learn More](#)

Parameter	Value
Product Name	PTBP2 Rabbit mAb
Remarks/Alias	nPTB; PTBLP; brPTB; PTBP2
Species	Human
Gene ID (Human)	58155
Gene ID	58155
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PTBP2 (Q9UKA9)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKA9
Protein Weight	57kDa
Shipping	Ice bag