

Anti-Ptcd1 Polyclonal Antibody For Wb, If/Icc, Elisa - 075127

Item Number: CM0021074

Introduction

Rabbit polyclonal antibody against human PTC1, suitable for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Designed for mitochondrial tRNA processing and regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PTCD1 Rabbit pAb
Remarks/Alias	PTCD1
Species	Human
GeneID	26024
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 435-537 of human PTC1 (NP_056360.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075127
Protein Weight	79kDa
Shipping	Ice bag