


Anti-Ptcd3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96Ey7

Item Number: CM0026838



Introduction

Rabbit polyclonal antibody against human PTCD3 (mS39). Validated for WB, IF/ICC, ELISA. Detects mitochondrial small ribosomal subunit protein mS39 in human, mouse, and rat samples. Suitable for mitochondrial research and disease diagnostics.

[Learn More](#)

Parameter	Value
Product Name	PTCD3 Rabbit pAb
Remarks/Alias	COXPD51; MRP-S39; PTCD3
Species	Human
GeneID (Human)	55037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 450-689 of human PTCD3 (NP_060422.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96EY7
Protein Weight	79kDa
Shipping	Ice bag