


Anti-Ptch1 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q13635

Item Number: CM0013757



Introduction

Rabbit polyclonal antibody against human PTCH1, a hedgehog receptor and tumor suppressor. Suitable for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for signal transduction and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PTCH1 Rabbit pAb
Remarks/Alias	PTC; BCNS; PTC1; PTCH; NBCCS; PTCH1
Species	Human
GeneID (Human)	5727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1214-1312 of human PTCH1 (NP_000255.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13635
Protein Weight	161kDa
Shipping	Ice bag