

Anti-Pten Monoclonal Antibody For Wb, Elisa - P60484

Item Number: CM0002456



Introduction

High-quality mouse monoclonal antibody targeting human PTEN (P60484), suitable for Western blot and ELISA. Cross-reacts with mouse and rat PTEN. Ideal for cancer research and PI3K-AKT pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PTEN Mouse mAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; PTEN
Species	Human
Gene ID (Human)	5728
Gene ID	5728
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 204-403 of human PTEN (P60484))
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P60484
Protein Weight	47 kDa
Shipping	Ice bag