

Anti-Pten Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P60484

Item Number: CM0002449



Introduction

KO-validated rabbit monoclonal antibody against human PTEN (P60484). Suitable for Western blot, immunoprecipitation, ELISA, and IF/ICC. Cross-reacts with mouse and rat. Ideal for cancer and signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PTEN Rabbit mAb
Remarks/Alias	BZS, DEC, CWS1, GLM2, MHAM, TEP1, MMAC1, PTEN1, 10q23del, PTENbeta, EN
Species	Human
Gene ID	5728 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-403 of human PTEN (P60484)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag