

Anti-Pten Monoclonal Antibody For Wb, Ihc-P, Elisa - P60484

Item Number: CM0004427



Introduction

Rabbit monoclonal antibody against human PTEN (phosphatase and tensin homolog) for WB, IHC-P, ELISA. Validated in human, mouse, rat. Targets C-terminal region (aa304-403). Ideal for cancer and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	PTEN Rabbit mAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; PTEN
Species	Human
Gene ID (Human)	5728
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 304-403 of human PTEN (P60484)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47 kDa
Shipping	Ice bag