

Anti-Pten Polyclonal Antibody For Ihc-P, If/Icc, Elisa - P60484

Item Number: CM0023151



Introduction

CamelBio offers rabbit polyclonal anti-PTEN antibody (P60484) validated for IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat PTEN, a key PI3K-AKT antagonist and tumor suppressor.

[Learn More](#)

Parameter	Value
Product Name	PTEN Rabbit pAb
Remarks / Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; PTEN
Species	Human
Gene ID (Human)	5728
Gene ID	5728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-304 of human PTEN (NP_000305.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag