

Anti-Pten Rabbit Monoclonal Antibody For Wb - P60484

Item Number: CM0005445

Introduction

KO-validated rabbit anti-PTEN monoclonal antibody for Western blot, IP, IF/ICC, IHC-P, and ELISA. Targets human PTEN (P60484) with cross-reactivity to mouse and rat. Ideal for oncology and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] PTEN Rabbit mAb
Remarks/Alias	BZS, DEC, CWS1, GLM2, MHAM, TEP1, MMAC1, PTEN1, 10q23del, PTENbeta, PTENGama
Species	Human
Gene ID (Human)	5728
Gene ID	5728
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 304-403 of human PTEN (P60484)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	20 kDa/47 kDa/65 kDa
Shipping	Ice bag