

Anti-Pten Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P60484

Item Number: CM0014344



Introduction

PTEN Rabbit Polyclonal Antibody (P60484) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets the tumor suppressor PTEN, a dual-specificity phosphatase regulating the PI3K/AKT pathway.

[Learn More](#)

Parameter	Value
Product Name	PTEN Rabbit pAb
Remarks/Alias	BZS; DEC; CWS1; GLM2; MHAM; TEP1; MMAC1; PTEN1; 10q23del; PTENbeta; PTEN
Target Species	Human
Gene ID (Human)	5728
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human PTEN (NP_000305.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60484
Protein Weight	47kDa
Shipping	Ice bag