

Anti-Ptgds Polyclonal Antibody For Wb, If/Icc, Elisa - P41222

Item Number: CM0021596



Introduction

Rabbit polyclonal antibody targeting human PTGDS (Lipocalin-type prostaglandin-D synthase). Validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for prostaglandin metabolism and CNS research.

[Learn More](#)

Parameter	Value
Product Name	PTGDS Rabbit pAb
Remarks/Alias	PDS; PGD2; PGDS; LPGDS; PGDS2; L-PGDS; PTGDS
Species	Human
GeneID (Human)	5730
GeneID	5730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-190 of human PTGDS (NP_000945.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41222
Protein Weight	21kDa
Shipping	Ice bag