

Anti-Ptges2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9H7Z7

Item Number: CM0028306

Introduction

Rabbit polyclonal antibody targeting human PTGES2 (mPGES-2), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for prostaglandin metabolism and cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	PTGES2 Rabbit pAb
Remarks/Alias	GBF1; GBF-1; PGES2; C9orf15; mPGES-2; PTGES2
Species	Human
Gene ID (Human)	80142
Gene ID	80142
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 88-377 of human PTGES2 (NP_079348.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H7Z7
Protein Weight	42kDa
Shipping	Ice bag