

Anti-Ptgfr Rabbit Polyclonal Antibody For Wb - P43088

Item Number: CM0015074



Introduction

Rabbit polyclonal antibody targeting human PTGFR (Prostaglandin F2-alpha receptor). Validated for WB and ELISA, with cross reactivity to mouse and rat. Ideal for GPCR signaling research.

[Learn More](#)

Parameter	Specification
Product Name	PTGFR Rabbit pAb
Remarks / Alias	FP; PTGFR
Species	Human
GeneID (Human)	5737
GeneID	5737
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P43088
Protein Weight	40kDa
Shipping	Ice bag