


Anti-Ptgir Polyclonal Antibody For Wb, Ihc-P, Elisa - P43119

Item Number: CM0019681



Introduction

PTGIR Rabbit pAb validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Detects prostacyclin receptor, a key GPCR target in cardiovascular and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	PTGIR Rabbit pAb
Remarks/Alias	IP; PRIPR; PTGIR
Species	Human
GeneID (Human)	5739
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 207-386 of human PTGIR (NP_000951.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43119
Protein Weight	41kDa
Shipping	Ice bag