

Anti-Ptgr3/Zadh2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N4Q0

Item Number: CM0020967



Introduction

High-performance rabbit polyclonal antibody against PTGR3/ZADH2, validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. 40 kDa. Suitable for prostaglandin reductase and adipocyte differentiation research.

[Learn More](#)

Parameter	Value
Product Name	PTGR3/ZADH2 Rabbit pAb
Remarks/Alias	PRG-3; ZADH2; PTGR-3; PTGR3/ZADH2
Species	Human
GeneID (human)	284273
GeneID	284273
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 33-371 of human PTGR3/ZADH2 (NP_787103.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4Q0
Protein Weight	40kDa
Shipping	Ice bag