


Anti-Pth1R Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q03431

Item Number: CM0023291



Introduction

Rabbit polyclonal antibody against human PTH1R (parathyroid hormone 1 receptor), suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying PTH/PTHrP signaling and bone/kidney receptor function.

[Learn More](#)

Parameter	Value
Product Name	PTH1R Rabbit pAb
Remarks/Alias	PFE, EKNS, PTHR, PTHR1, PTH1R
Species	Human
Gene ID (Human)	5745
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-188 of human PTH1R (NP_000307.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03431
Protein Weight	66kDa
Shipping	Ice bag