


Anti-Pthlh Polyclonal Antibody For Wb, Ihc-P, Elisa - P12272

Item Number: CM0029100



Introduction

Rabbit polyclonal antibody targeting human PTHLH. Validated for WB, IHC-P, ELISA with reactivity in human, mouse, rat. Ideal for bone development, calcium transport, and cancer cell migration research.

[Learn More](#)

Parameter	Value
Product Name	PTHLH Rabbit pAb
Remarks/Alias	HHM, PLP, BDE2, PTHR, PTHRP, PTHLH
Species	Human
Gene ID (Human)	5744
Gene ID	5744
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-175 of human PTHLH (NP_945316.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12272
Protein Weight	20kDa
Shipping	Ice bag