


Anti-Pthlh Rabbit Polyclonal Antibody For Wb, Elisa - P12272

Item Number: CM0021257



Introduction

Polyclonal rabbit antibody against human PTHLH, validated for WB and ELISA. Reacts with human, mouse, rat. Ideal for studying parathyroid hormone-related protein in development and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	PTHLH Rabbit pAb
Remarks/Alias	HHM; PLP; BDE2; PTHR; PTHRP; PTHLH
Species	Human
Gene ID (Human)	5744
Gene ID	5744
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-175 of human PTHLH (NP_945315.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12272
Protein Weight	20 kDa
Shipping	Ice bag