


Anti-Ptk2B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14289

Item Number: CM0020322



Introduction

Explore PTK2B (FAK2/PYK2) rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. For research on focal adhesion kinase signaling, immune cell migration, and bone remodeling.

[Learn More](#)

Parameter	Value
Product Name	PTK2B Rabbit pAb
Remarks/Alias	PKB; PTK; CAKB; FAK2; PYK2; CADTK; FADK2; RAFTK; PTK2B
Species	Human
GeneID	2185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 780-1009 of human PTK2B (NP_004094.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14289
Protein Weight	116kDa
Shipping	Ice bag