

Anti-Ptk6 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13882

Item Number: CM0030340



Introduction

Rabbit polyclonal antibody against human PTK6 (BRK). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for epithelial biology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PTK6 Rabbit pAb
Remarks/Alias	BRK; PTK6
Species	Human
Gene ID (Human)	5753
Gene ID	5753
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 222-451 of human PTK6 (NP_005966.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13882
Protein Weight	52kDa
Shipping	Ice bag