


Anti-Ptk7 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13308

Item Number: CM0030150



Introduction

Rabbit polyclonal antibody targeting human PTK7/CCK4, validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for Wnt signaling and cell adhesion research.

[Learn More](#)

Parameter	Specification
Product Name	PTK7 Rabbit pAb
Remarks/Alias	CCK4; CCK-4; PTK7
Species	Human
GeneID (Human)	5754
GeneID	5754
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 811-1070 of human PTK7 (NP_002812.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13308
Protein Weight	118kDa
Shipping	Ice bag