


Anti-Ptp1B Monoclonal Antibody For Wb, Ihc-P, Elisa - P18031

Item Number: CM0007166



Introduction

High-quality rabbit monoclonal antibody against human PTP1B, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for tyrosine phosphatase signaling and ER stress pathway research.

[Learn More](#)

Parameter	Value
Product Name	PTP1B Rabbit mAb
Remarks/Alias	PTP1B
Species	Human
GeneID (human)	5770
GeneID	5770
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 281-435 of human PTP1B (NP_002818.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P18031
Protein Weight	50kDa
Shipping	Ice bag