


Anti-Ptp1B (Ptpn1) Polyclonal Antibody For Wb, If/Icc, Elisa - P18031

Item Number: CM0018856



Introduction

Rabbit polyclonal antibody against human PTPN1 (PTP1B) validated for WB, IF/ICC, and ELISA, cross-reacting with mouse and rat. Useful for studying tyrosine phosphatase roles in ER stress response, signal transduction, and cell migration. Ships on ice.

[Learn More](#)

Parameter	Value
Product Name	PTPN1 Rabbit pAb
Remarks/Alias	PTP1B; PTPN1
Species	Human
GeneID (Human)	5770
GeneID	5770
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 281-435 of human PTPN1 (NP_002818.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18031
Protein Weight	50kDa
Shipping	Ice bag