

Anti-Ptp4A3 Polyclonal Antibody For Wb, Ihc-P, Elisa - 075365

Item Number: CM0022633



Introduction

Rabbit polyclonal antibody targeting human PTP4A3 (PRL-3), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Ideal for cancer metastasis, cell proliferation, and cardiac hypertrophy research.

[Learn More](#)

Parameter	Value
Product Name	PTP4A3 Rabbit pAb
Remarks/Alias	PRL3; PRL-3; PRL-R; PTP4A3
Species	Human
GeneID	11156
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-148 of human PTP4A3 (NP_009010.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	075365
Protein Weight	20kDa
Shipping	Ice bag