


Anti-Ptpip51/Rmdn3 Monoclonal Antibody For Wb, Elisa - Q96Tc7

Item Number: CM0007346



Introduction

Rabbit monoclonal antibody against human PTPIP51/RMDN3, validated for WB and ELISA, with cross-reactivity to rat. Suitable for studying calcium homeostasis, apoptosis, and mitochondrial dynamics.

[Learn More](#)

Parameter	Value
Product Name	PTPIP51/RMDN3 Rabbit mAb
Remarks/Alias	RMD3; RMD-3; FAM82C; FAM82A2; ptpip51; PTPIP51/RMDN3
Species	Human
GeneID (Human)	55177
GeneID	55177
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human PTPIP51/RMDN3 (NP_060615.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q96TC7
Protein Weight	52kDa
Shipping	Ice bag