


Anti-Ptpip51/Rmdn3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Tc7

Item Number: CM0027730



Introduction

Rabbit polyclonal antibody against PTPIP51/RMDN3, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Useful for studying calcium homeostasis, apoptosis, and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	PTPIP51/RMDN3 Rabbit pAb
Remarks/Alias	RMD3; RMD-3; FAM82C; FAM82A2; ptpip51; PTPIP51/RMDN3
Species	Human
Gene ID	55177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 171-470 of human PTPIP51/RMDN3 (NP_060615.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96TC7
Protein Weight	52kDa
Shipping	Ice bag